

IHE Integration Statement			
Vendor	Product Name	Version	Date
Toshiba Medical Systems Corporation	Aquilion™ Aquilion ONE™ Aquilion™ PRIME Alexion™ TSX-101A,TSX-301A, TSX-301B,TSX-301C, TSX-302A,TSX-032A,TSX-201A	V4.93 or Later	1 Oct.2013
This product implements all transactions required in the IHE Technical Framework (Rev 12.0)and IHE Cardiology Technical Framework (Rev 5.0) to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol Setting Modality Group Case PPS Exception Management Billing and Material Management Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Charge Posting	Acquisition Modality	PPS Exception Management Modality Group Case	
Presentation of Grouped Procedures	Acquisition Modality	No options defined	
Evidence Documents	Acquisition Modality	TID 3900 Cardiovascular Analysis Report Calcium Score Evidence	
Portable Data for Imaging	Portable Media Creator	None	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Cardiology			
Evidence Documents	Acquisition Modality	TID 3900 Cardiovascular Analysis Report Calcium Score Evidence	
Internet address for vendor's IHE information:			
Japan: http://www.toshiba-medical.co.jp/tmd/products/ihe http://www.toshiba-medical.co.jp/tmd/english/ihe (International page)			
North America: http://medical.toshiba.com/services/enterprise-integration/index.php			
Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements			

Links to Standards Conformance Statements for the Implementation		
HL7	None	
DICOM	Japan: http://www.toshiba-medical.co.jp/tmd/products/dicom http://www.toshiba-medical.co.jp/tmd/english/dicom (International page) North America: http://medical.toshiba.com/services/enterprise-integration/index.php Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements	
Links to general information on IHE		
In North America: http://www.ihe.net/	In Europe: http://www.ihe-europe.net/	In Japan: http://www.ihe-j.org/

“Aquilion”, “Alexion” and “Aquilion ONE” is a trademark of Toshiba Medical Systems Corporation.