

IHE Integration Statement			
Vendor	Product Name	Version	Date
Toshiba Medical Systems Corporation	Aplio™ XG SSA-790A Aplio™ MX SSA-780A Xario™ XG SSA-680A Xario™ SSA-660A	V5.20 V2.20 V4.20 V9.20	March 2012
This product implements all transactions required in the IHE Radiology Technical Framework (Rev 10.0), IHE Cardiology Technical Framework (Rev 2.1), and IHE IT Infrastructure Technical Framework (Rev 8.0) to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Evidence Documents	Acquisition Modality	No options defined TID 5000 OB-GYN Ultrasound Procedure Report	
Teaching file and Clinical trial Export	Export Selector	Additional Teaching File Information Delay for Reason	
Cardiology			
Echocardiography Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo	
Evidence Documents	Acquisition Modality	Echo Evidence TID 5100 Vascular Ultrasound Report TID 5200 Echocardiography Procedure Report	
Stress Testing Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo	
IT Infrastructure			
Consistent Time	Time Client	None	

<u>Internet address for vendor's IHE information:</u>		
Japan: http://www.toshiba-medical.co.jp/tmd/products/ihe http://www.toshiba-medical.co.jp/tmd/english/ihe (International page)		
North America: http://medical.toshiba.com/services/enterprise-integration/index.php		
Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/Ultrasound/Compliance-Statements		
Links to Standards Conformance Statements for the Implementation		
HL7	None	
DICOM	Japan: http://www.toshiba-medical.co.jp/tmd/products/dicom http://www.toshiba-medical.co.jp/tmd/english/dicom (International page) North America: http://medical.toshiba.com/services/enterprise-integration/index.php Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/Ultrasound/Compliance-Statements	
Links to general information on IHE		
In North America: http://www.ihe.net/	In Europe: http://www.ihe-europe.net/	In Japan: http://www.ihe-j.org/

"Aplio" and "Xario" are trademarks of Toshiba Medical Systems Corporation.