

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Toshiba Medical Systems Corporation	Aplio™ 500 TUS-A500 Aplio™ 400 TUS-A400 Aplio™ 300 TUS-A300	V3.00	February 2014
This product implements all transactions required in the IHE Radiology Technical Framework (Rev 12.0), IHE Cardiology Technical Framework (Rev 5.0) and IHE IT Infrastructure Technical Framework (Rev 10.1) to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b>			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Assisted Acquisition Protocol Setting Modality Group Case	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Evidence Documents	Acquisition Modality	No options defined TID 5000 OB-GYN Ultrasound Procedure Report	
<b>Cardiology</b>			
Echocardiography Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo	
Evidence Documents	Acquisition Modality	Echo Evidence TID 5100 Vascular Ultrasound Report TID 5200 Echocardiography Procedure Report	
Stress Testing Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo	
<b>IT Infrastructure</b>			
Consistent Time	Time Client	None	

<b>Internet address for vendor's IHE information:</b>		
Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/ihe">http://www.toshiba-medical.co.jp/tmd/products/ihe</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/ihe">http://www.toshiba-medical.co.jp/tmd/english/ihe</a> (International page)		
North America: <a href="http://medical.toshiba.com/services/enterprise-integration/index.php">http://medical.toshiba.com/services/enterprise-integration/index.php</a>		
Europe: <a href="http://www.toshiba-medical.eu/en/Our-Product-Range/Ultrasound/Compliance-Statements">http://www.toshiba-medical.eu/en/Our-Product-Range/Ultrasound/Compliance-Statements</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>HL7</b>	None	
<b>DICOM</b>	Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/dicom">http://www.toshiba-medical.co.jp/tmd/products/dicom</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/dicom">http://www.toshiba-medical.co.jp/tmd/english/dicom</a> (International page) North America: <a href="http://medical.toshiba.com/services/enterprise-integration/index.php">http://medical.toshiba.com/services/enterprise-integration/index.php</a> Europe: <a href="http://www.toshiba-medical.eu/en/Our-Product-Range/Ultrasound/Compliance-Statements">http://www.toshiba-medical.eu/en/Our-Product-Range/Ultrasound/Compliance-Statements</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Europe: <a href="http://www.ihe-europe.net/">http://www.ihe-europe.net/</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>

"Aplio" is a trademark of Toshiba Medical Systems Corporation.