

IHE Integration Statement			
Vendor	Product Name	Version	Date
Toshiba Medical Systems Corporation	Infinix™ i-series / DFP-8000	V4.25	25 Jan. 2011
This product implements all transactions required in the IHE Technical Framework (Rev 9.0), IHE Cardiology Technical Framework (Rev 2.1) and IHE IT Infrastructure Technical Framework (Rev 3.0) to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Cardiology			
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query	
IT Infrastructure			
Consistent Time	Time Client	None	
<u>Internet address for vendor's IHE information:</u>			
Japan: http://www.toshiba-medical.co.jp/tmd/products/ihe http://www.toshiba-medical.co.jp/tmd/english/ihe (International page)			
North America: http://medical.toshiba.com/service/ihe/ihe-xr.php			
Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/X-ray/Compliance-Statements			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: http://www.toshiba-medical.co.jp/tmd/products/dicom http://www.toshiba-medical.co.jp/tmd/english/dicom (International page) North America: http://medical.toshiba.com/service/dicom/dicom-xr.php Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/X-ray/Compliance-Statements		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

"Infinix" is a trademark of Toshiba Medical Systems Corporation.