

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Toshiba Medical Systems Corporation	DI Station™ ADR-1000A	V3.1	3 Feb. 2006
This product implements all transactions required in the IHE Technical Framework (Rev 6.0) to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b>			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
<b>Internet address for vendor's IHE information:</b>			
Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/ihe/">http://www.toshiba-medical.co.jp/tmd/products/ihe/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/ihe/">http://www.toshiba-medical.co.jp/tmd/english/ihe/</a> (International page)			
North America: <a href="http://medical.toshiba.com/services/iheintegrationstatements.htm">http://medical.toshiba.com/services/iheintegrationstatements.htm</a>			
Europe: <a href="http://www.toshiba-europe.com/medical/medical.asp?PageID=153">http://www.toshiba-europe.com/medical/medical.asp?PageID=153</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/dicom/">http://www.toshiba-medical.co.jp/tmd/products/dicom/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/dicom/">http://www.toshiba-medical.co.jp/tmd/english/dicom/</a> (International page)		
	North America: <a href="http://medical.toshiba.com/services/dicomconformancstatements.htm">http://medical.toshiba.com/services/dicomconformancstatements.htm</a>		
	Europe: <a href="http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532">http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Europe: <a href="http://www.ihe-europe.org/">http://www.ihe-europe.org/</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>	

“DI Station“ is a trademark of Toshiba Medical Systems Corporation.